Aqua-Pure

Spec #

Quantity

Models:

AP Easy CS-FF / AP Easy Cyst-FF

Drinking Water Systems

Applications:

- Chlorine Taste & Odor
- Sediment*
- Parasitic Protozoan Cysts (AP Easy Cyst-FF)

At a Glance:

- · Easy to select the right product
- · Easy to install, no special tools required
- Easy to maintain, filter change is a twist





AP Easy CS-FF

Product Benefits:

- Designed to hook up to your existing kitchen or bath faucet for full flow filtration, providing a clean look on your counter
- Both systems reduce CTO and Sediment. The AP Easy Cyst-FF removes 99.95% of filterable Parasitic Protozoan Cysts from your cold water supply right at the sink
- Everything included no additional parts necessary for most installations

AP Easy Cyst-FF


Product Specifications:


Model No. (Part No.)	Inlet/ Outlet Size	Dimensions				Мах.	Max.	Std.
		Height	Diameter	Flow Rate	Capacity	Water Temp.	Water Pressure	Replacement Ctg. No.
AP Easy CS-FF (56092-22)	3/8" Push	12" (30.5 cm)	3-3/16" (8.1 cm)	2.0 gpm (7.8 lpm)	2,000 gallons (7,571 liters)	100°F (37.8°C)	125 psi (862 kPa)	C-CS-FF 56104-27
AP Easy Cyst-FF (56092-23)		15 5/32" (38.5 cm)	3-3/16" (8.1 cm)	1.5 gpm (5.7 lpm)				C-Cyst-FF 56104-28

* Claims not tested or certified by NSF









Spec #

Quantity

Model Number

Part Number

Physical Specifications:

- 3/8" Push-In fittings
- · System is a single vessel design
- System is manufactured from FDA CFR-21 compliant materials
- System has a maximum operating pressure of 125 psi (862 kPa) and operating temperature of 100°F (37.8°C)
- · System includes attached mounting bracket
- System to be installed to primary kitchen or bath cold water faucet

WARNING: Do not use with water that is microbiologically unsafe or of unknown quality without adequate disinfection before or after the system. The contaminants or other substances removed or reduced by this water filter are not necessarily in all users water. Systems certified for cyst reduction may be used on disinfected water that may contain filterable cysts.

EPA Establishment Number 070595-CT-001

CUNO Incorporated recommends regularly scheduled maintenance and replacement of the filter cartridge(s) in order for the product to perform as advertised/sold. Change the filter cartridge at least every 6 months. Improper installation and maintenance may result in property damage due to water leakage.

Limited Warranty:

CUNO Incorporated warrants that this product is free from defects in materials and workmanship. This limited warranty, together with any and all warranties implied by law, shall be limited to a period of one (1) year from the date of original purchase. This limited warranty does not apply to failures that result from abuse, misuse, alteration or damage not caused by CUNO or failure to properly comply with these installation or cartridge change-out instructions. This limited warranty is subject to exclusions and limitations. Please refer to the Warranty Card provided with the product for details.

Distributed by:



Visit www.nsf.org for the claims associated with products that are NSF listed.





CUNO Incorporated
400 Research Parkway
Meriden, CT 06450 U.S.A.
Toll Free: 1.800.222.7880
Worldwide: 203.237.5541
Fax: 203.238.8701
www.aguapure.com = www.cuno.com