

Quatromyicide V

Disinfectant, No Rinse Sanitizer, and Deodorizer



Applications:

Disinfection:

To disinfect pre-cleaned, inanimate, hard surfaces apply this product with mop, cloth, sponge or coarse (trigger) spray so as to wet all surfaces thoroughly. Allow to remain wet for 10 minutes, then remove excess liquid. Prepare a fresh solution for each use.

Sanitizing:

Remove all gross food particles and soil from areas which are to be sanitized with good detergent, pre-flush, pre-soak or pre-scrape treatment. Rinse with a potable water rinse. To sanitize pre-cleaned and potable water-rinsed, non-porous, food contact surfaces, prepare a 200 ppm active quaternary solution by adding 8 mL of this product to 4 L of water.

Benefits:

- Disinfects, sanitizes, and deodorizes
- Great solution for multi-touch surfaces responsible for cross-contamination
- Deodorizes those areas which are generally hard to keep fresh smelling
- Ideal for the food industry, schools, hospitals and delicate personnel settings such as daycares or senior care facilities
- No rinse food contact sanitizer when used at 200 ppm

Specifications:

| | |
|--------------|----------------|
| Colour | Colourless |
| DIN | 02507064 |
| Contact Time | 10 minutes |
| Odour | Mild |
| pH | 7.0 - 8.0 |
| Viscosity | Not determined |
| Solvents | No |
| Phosphates | No |
| WHMIS | Exempt |
| ECOLOGO® | No |
| Food Safe | Yes |
| Kosher | Yes |
| Shelf Life | Expiry date |

Codes:

- **Product Code:** 52898
- **UPC:** 0 56069 52898 6

Dilution Rates:

* Tip colour proportioning systems at 3.5 GPM

General Disinfection /
Hospital Disinfection /
Deodorization / Mildewstat

22 mL / 4 L (1:160)



Sanitization

8 mL / 4 L (1:500)

Available Format:



4 L

Packaging:

Bottle: 5.8" x 5.8" x 12.25"
(14.7 x 14.7 x 31 cm)
10 lbs (4.5 kg)

Box: 11.625" x 11.625" x 12.875"
40 lbs (18 kg)
4 per case

